

cura¹
Fall Prevention Products



We have a range of monitors to alert when a wireless device, such as a pendant or chair pad, is activated in your communal areas.

It is like having an extra pair of eyes on all your falls risk residents / patients.



SOS Buttons



Chair Pads



Motion Sensors



Wearable Pendants



Communal Area Wireless Monitoring

Be alerted when a falls risk resident / patient gets up from their chair or wheelchair in your communal areas.

Also allow them to push a button to notify that assistance is required.

Our 6 channel monitor can also plug into your nurse call system via a connection lead.

*please note: you may need to organise for your nurse call provider to change the call point in the communal area to one with a socket.



Type	Picture	Best Features	Nurse Call Connection	# Pairable Devices	Portable	Range*
6 Channel Monitor		Can monitor chair pads and wireless pendants. Will emit a volume-controlled alarm and/or go through the nurse call system (contact us for more information) via a nurse call connection lead.	Yes	6	Yes	30-50m
ActiveCare LCD Pager		This pager can work as a standalone alert if people get up from chairs or press pendants. Can also work in conjunction with the 6 Channel Monitor.	No	9	Yes	30-50m

*Listed range is line of sight and subject to environmental conditions. We recommend testing range before use.



Qube Industries
Unit 1, 11 Ring Street, Neerabup, WA 6031
admin@qubeindustries.com.au 1300 577 757

HOW TO ORDER:












Monitors: Select the monitor that best suits your facilities patients/residents requirements.

Picture	Code	Description	More information
	2700	6 Channel Wireless Monitor	Adjustable alarm tone & volume. Alarm memory & automatic reset. Low battery & lost pad indicator. 3 x AA internal battery or external 6Vdc plug pack.
	3501	ActiveCare LCD Pager	Large 7 Character alpha numeric display. Red reminder LED. Selectable OFF/Vibrate/Vibrate & Tone. Selectable alarm cancel. Optional wall mount bracket. Includes belt clip and batteries.

Connection Lead: If using a 2700 6 Channel Monitor and you would like to connect this to a call point you will require a nurse call connection lead. *Please note: if the call point in your communal area does not have a socket you will need to organise with your nurse call provider to have this changed.

Picture	Code	Description	More information
	2460	Nurse Call Connection Lead	6.35mm stereo to 6.35mm stereo Plug (TS) – 2m
	2461	Nurse Call Connection Lead	6.35mm mono to 6.35mm stereo (TS) – 2m
	2643	Nurse Call Connection Lead	6.35mm mono to 6.35mm mono (TS) – 2m – suits most systems
	2646	Nurse Call Connection Lead	6.35mm mono to 6.35mm mono (TS) – 3m – suits most systems
	2493	Nurse Call Connection Lead	6.35mm mono to 3.5mm mono (TS) – 2m

Wireless Programmable Transmitters: Choose the devices to best suit requirements. Add devices at any stage. All pads require a 2719 or 2720 transmitter.

Picture	Code	Description	Best Feature	Transmitter Required
	2630	Wireless Chair Pad	Chair pad can be placed anywhere in the room	2719
	2691	Universal Chair Pad	18mth battery life/warranty	2720
	2719	Replacement Transmitter	12-15mth battery life	-
	2720	Replacement Universal Transmitter	18mth battery life/warranty	-
	2721	PIR Motion Sensor	Ideal as an out of/into room sensor	-
	3503	Door Reed Switch	Ideal as an out of/into room sensor	-
	2588	Pendant with Lanyard	Economical & lightweight	-
	2595	Pendant on Wristband	Low battery warning LED	-
	2599	SOS Button	Large button, easy to see and use	-
	3502	Handset with Reset	Easy to use wireless handset	-
	2558	Waterproof Call Point with Reset	Pull string & push button activation	-