

**cura<sup>1</sup>**  
**Fall Prevention Products**



**The cura<sup>1</sup> Pull String Monitor is simple but effective.**

**Pull string activates the alarm.**

**Simply clip the pull-string to the person's clothing to raise the alarm when they move away.**



Walker Monitoring



Out Of Chair Monitoring



Wheelchair Monitoring

**2617 Pull String Monitor**



**Features:**

- \*Simple operation and maintenance.
- \*Affordable price.
- \*Pull string activates the alarm.
- \*Auto re-arm and alarm memory.
- \*Battery operated (battery included).
- \*Adjustable alarm volume.
- \*Red LED alarm indicator light.
- \*Low battery detection and indication.

**Range of uses:**

- \*Connect the Pull String Monitor to a walker to alert when a person moves away from their walker.
- \*Connect the Pull String Monitor to a wheelchair to alert when a person gets up.
- \*Connect the Pull String Monitor to a chair to alert when a person stands up.

**Qube Industries**  
**Unit 1, 11 Ring Street, Neerabup, WA 6031**  
**[admin@qubeindustries.com.au](mailto:admin@qubeindustries.com.au) 1300 577 757**