



**Our ActiveCare LCD Pager can be standalone or work in conjunction with the ActiveCare console or 6 Channel Monitor.**

**It can be programmed with up to 9 devices.**

**Ideal for facility and home use.**

### **3501 ActiveCare LCD Pager**

#### **Features:**

- \*Large easy-to-read 7-digit LCD screen displays device type and ID.
- \* Can be programmed to monitor up to 9 wireless devices.
- \*Up to 5 simultaneous alarms can be stored in the pager's "Call Review" memory.
- \*Vibration and tone selection for hard of hearing.
- \*Out of range check.
- \*Wireless reset capability.
- \*Direct signalling from wireless devices for use in home care environments.
- \*Easily cleanable buttons for a hospital setting.
- \*Typical operating range: 30 – 50m (dependent on environmental conditions).
- \*Power requirements: 2 x AA batteries (included).
- \*Dimensions: 83 x 55 x 26mm.
- \*Housing material: Durable ABS plastic.



#### **Range of uses:**

- \*The ActiveCare Pager will alert when a falls risk person is at risk or requires attention depending on the wireless devices programmed. It can be used for residential facilities or home care use.
- \*The pager can be used as a standalone alert in communal areas of residential facilities for people getting up from chairs or pressing their pendants.
- \*The pager can be used to extend the range of the 6 Channel Monitor (2700) or the ActiveCare Console (3500).

**Qube Industries**  
**Unit 1, 11 Ring Street, Neerabup, WA 6031**  
**[admin@qubeindustries.com.au](mailto:admin@qubeindustries.com.au) 1300 577 757**

## HOW TO ORDER:

Start with the 3501 ActiveCare LCD Pager and then choose up to 9 wireless programmable transmitters from the range listed below.

Picture	Code	Description
	3501	ActiveCare LCD Pager

**Wireless Programmable Transmitters:** Choose the devices to best suit requirements. You can choose up to 9 devices. Add devices at any stage. All pads require a 2719 or 2720 transmitter.

Picture	Code	Description	Best Feature	Transmitter Required
	2630	Wireless Chair Pad	Chair pad can be placed anywhere in the room	2719
	2691	Universal Chair Pad	18mth battery life/warranty	2720
	2629	Wireless Above Mattress Bed Pad	Can detect people 40kg & under	2719
	3515	Wireless Under Mattress Bed Pad	People won't even know it is there	2719
	2670	Wireless PreAlert Bed Pad	Alerts when they move to the edge of the bed	-
	2690	Universal Above Mattress Bed Pad	18mth battery life/warranty	2720
	2624	Wireless Floor Sensor Mat	Can be moved around the room as required	2719
	2671	Universal Single Fold Floor Sensor Mat	18mth battery life/warranty	2720
	2672	Universal Dual Fold Floor Sensor Mat	18mth battery life/warranty	2720
	2692	Universal Floor Sensor Mat	18mth battery life/warranty	2720
	2719	Replacement Transmitter	12-15mth battery life	-
	2720	Replacement Universal Transmitter	18mth battery life/warranty	-
	2721	PIR Motion Sensor	Ideal as an out of/into room sensor	-
	3503	Door Reed Switch	Ideal as an out of/into room sensor	-
	2588	Pendant with Lanyard	Economical & lightweight	-
	2595	Pendant on Wristband	Low battery warning LED	-
	2599	SOS Button	Large button, easy to see and use	-
	3502	Handset with Reset	Easy to use wireless handset	-
	2558	Waterproof Call Point with Reset	Pull string & push button activation	-